


POCKET MICROSCOPE

NATURE DISCOVERY KIT



What's inside your experiment kit:



>>> KIT CONTENTS

Checklist: Find – Inspect – Check off

✓	No.	Description	Qty.	Item No.
<input type="radio"/>	1	Left housing	1	718 024
<input type="radio"/>	2	Right housing	1	718 026
<input type="radio"/>	3	Housing cover portion	1	718 027
<input type="radio"/>	4	Microscope stand	1	718 176
<input type="radio"/>	5	Transparent bottom housing attachment	1	718 029
<input type="radio"/>	6	Zoom adjustment wheel	1	718 030
<input type="radio"/>	7	Focus adjustment wheel	1	718 031
<input type="radio"/>	8	Eyepiece ring	1	718 033
<input type="radio"/>	9	Battery compartment cover	1	718 038
<input type="radio"/>	10	Eye cup	1	718 177
<input type="radio"/>	11	Eyepiece sleeve ring	1	718 178
<input type="radio"/>	12	Jacket sleeve	1	718 180
<input type="radio"/>	13	Middle sleeve with relay lens	1	718 181
<input type="radio"/>	14	Inner sleeve with objective lens	1	718 182
<input type="radio"/>	15	Button cell battery 1.5-volt, type LR 44	3	718 183
<input type="radio"/>	16	Prepared slide: <ul style="list-style-type: none"> · Fish scale · Snakeskin scale · Bird's feather 	1	718 043
<input type="radio"/>	17	Slide for your own specimens	1	718 044

✓	No.	Description	Qty.	Item No.
<input type="radio"/>	18	Illumination knob with large spring	1	718 032
<input type="radio"/>	19	Battery compartment closure knob with small spring	1	718 184
<input type="radio"/>	20	Eyepiece lens	1	718 185
<input type="radio"/>	21	Carrying strap	1	718 186
<input type="radio"/>	22	Electronic unit with switch, LED bracket, and wires	1	718 028
<input type="radio"/>	23	Screws for housing pieces	2	718 188
<input type="radio"/>	24	Screw for bottom housing attachment	1	718 187
<input type="radio"/>	25	Holding pin for carrying strap	1	718 189

Parts that are not included in the kit are indicated in *italic script* under the “YOU WILL NEED” heading in the experiments.

You will also need:

Small Phillips-head screwdriver, pen, tape, scissors (optional), interesting objects from around the house or from nature to study under the microscope

>>> TABLE OF CONTENTS

Kit Contents 1
Safety Information/A Word to Parents and Adults 2
Table of Contents 3

Your Pocket Microscope 4
 Assemble the microscope, insert the batteries, and off you go!



Assembly 5
First look 12

Exploring the World with your Microscope 14
 Here's where you will find suggestions for all the things you can study under your microscope.

Studying plants 15
Exciting bugs 16
Crime scene clues 17


TIP!

You will find additional information under "Check It Out" on pages 11, 18, 19, and 20.



ASSEMBLY


14




Insert the batteries. First stack the three button cells on top of one another, each with its positive end down. Then insert the stack into the battery compartment, again keeping the positive end down.

15

Mount the cover. First hook it onto the closure knob and then latch it in on the left side.



16



Done!

CHECK IT OUT




Light path

In putting together your microscope, you assembled the various lenses and their sleeves into a so-called “optical path.” The **eyepiece lens** sits fixed at the upper end. This is what your eye looks through.

The **relay lens** can be moved up and down with the zoom adjustment ring. It lets you switch between a low level of magnification (20 times) and high magnification (40 times).

The focus adjustment ring moves the **objective lens**. That adjusts the distance from the object in a way that lets you focus the image.



Circuit

When you installed the wires and electronic components, you built a circuit. From the positive terminal of the battery compartment, the red wire goes to the **switch** and then to the illumination unit, where the LED bulb is installed. From the bulb the black wire leads to the negative terminal of the battery compartment. When batteries are inserted and the switch is turned on, the circuit is closed. The current flows from the batteries through the **light**, and the LED lights up. If the switch is pushed again, the circuit is interrupted. Then the light goes out.



DID YOU KNOW?

An instrument that has at least a 40-fold magnification is called a microscope. Lower levels of magnification are often called “macroscopic.” Both words are based on Greek roots. “Mikros” means small, and “makros” means large. In either case, what it means is that small things can be seen as if they were noticeably larger.



EXPERIMENT 2



YOU WILL NEED

- › **Pocket microscope**
- › *Dead insect (from a window sill, for example)*

HERE'S HOW

Use your microscope without its stand. Place a dead insect on the table or another smooth surface and then position the microscope over it.

TIP!

There are lots of places where you can find “bugs” such as insects, spiders, or roly-polies. Take a look in the corners of the basement, under boards or rocks in the garden, or in a lampshade inside a room. Be sure to have your parents help you.

There are plenty of dead bugs that you will be able to find. Do not kill any bugs for your experiments!

BUGS UP CLOSE!

Insects and other bugs possess an **exoskeleton** made of chitin. On the surface, you will often find hairs, scales, or other structures. The eyes of insects consist of many individual eyes and are called **compound eyes**. On their feet, you can often see claws or suction disks, and the wings will have lots of fine veins running through them.

